**Absolute Kontraindikationen** Absolute contraindications Contraindicaciones absolutas:

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- Acute phase of myocardial infarction as well as severe cardiac insufficiency
- Acute leukemia and agranulocytosis
- Use in the case of bone fractures (exception: partially impacted wisdom teeth due to the risk of fracture gap infection)
- Dental ankylosis
- Radiation and chemotherapy

NOTE: Related contraindications and further instructions that are deemed set-specific and generally valid can be found at the download area of our website at "Application and Safety Instructions for Bone Management® Systems".

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## Inhalt Benex®-Control Professional

Content Contenido

### \* Benex®-Control Basic

|                     | Diamantierte Spiralbohrer<br>Diamond-coated twist drills<br>Taladros espirales diamantados |      | Extraktionsschrauben<br>Extraction screws<br>Tornillos de extracción |      |      |      | Schraubendreher<br>Extraction screw key<br>Atornilladores |       | Pilotbohrer <sup>3</sup><br>Pilot burs <sup>3</sup><br>Fresas piloto <sup>3</sup> |       |
|---------------------|--|------|--|------|------|------|---|-------|---|-------|
| Fig. BE001230...    | 030  | 035  | 060 *  | 070  | 080  | 085  | 045   | 055 * | A2001*  | B2001 |
| Shank <sup>1</sup>  | 206  | 206  |  |      |      |      |   |       | 205   | 205   |
| Size <sup>2</sup>   |  |      |  |      |      |      |   |       | 013   | 016   |
| Length mm           | 17.0   | 17.0 | 10.0   | 16.0 | 10.0 | 16.0 |   |       | 12.0  | 12.0  |
| □ Internal diameter | 1.3  | 1.8  | 1.3  | 1.3  | 1.8  | 1.8  |   |       |   |       |
| □ External diameter | 1.8  | 2.3  | 1.6  | 1.6  | 2.1  | 2.1  |   |       | 1.3   | 1.6   |

<sup>1</sup> 204=RA, 206=RA XL, 314=FG, 316=FG XL, 317=30 mm

<sup>2</sup>Largest working part diameter

<sup>3</sup>Alternative to diamond-coated twist drills

□ Internal diameter □ External diameter

|                     | Kronentrenner<br>Crown cutters<br>Fresas para separar coronas<br>(je 2 St. / 2 pcs. of each / 2 pcs. c/u) |        | Osteotomieinstrumente<br>Osteotomy instruments<br>Fresas para osteotomía |        |      |      | Schutzkappe<br>Protection cap<br>Casquillo protector |  |
|---------------------|---|--------|--|--------|------|------|--|--|
| Fig. BE001230...    | HM 31C  | HM 33C | HM331L   | HM341L | 859L | 859  | BE00123 SK60   |  |
| Shank <sup>1</sup>  | 204   | 314    | 316  | 317    | 204  | 204  |  |  |
| Size <sup>2</sup>   | 012   | 012    | 010  | 012    | 010  | 018  |  |  |
| Length mm           | 4.1   | 4.1    | 5.5  | 6.0    | 12.0 | 10.0 | 8.0  |  |
| □ Internal diameter |   |        |  |        |      |      | 5.0  |  |
| □ External diameter | 1.2   | 1.2    |  |        |      |      | 5.5  |  |

<sup>1</sup> 204=RA, 206=RA XL, 314=FG, 316=FG XL, 317=30 mm

<sup>2</sup>Largest working part diameter

□ Internal diameter □ External diameter

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<sup>3</sup>Teflon® is a registered trademark of DuPont, USA

|                  | Extraktor <sup>4</sup><br>Extractor <sup>4</sup><br>Extractor <sup>4</sup> | Zugseile <sup>4</sup><br>Traction strings <sup>4</sup><br>Cuerdas de tracción <sup>4</sup> | Quadrantenstütze <sup>4</sup><br>Quadrant Support <sup>4</sup><br>Soporte cuadrante <sup>4</sup> |
|------------------|--|--|--|
| Fig. BE001230... | ...010 *   | ...020 *   | ...090 *   |

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- Safe and controlled root extraction
- Maximum protection of hard bone and soft tissue
- Preservation of the vestibular bone lamella by means of axial force distribution
- Drilling also possible if only the smallest amount of root remains
- Reaming possible independent of the anatomical root canal
- Prevention of an osteotomy
- High level of patient acceptance
- Perfect alveolar management within the scope of an immediate implantation.

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Bone Management® is a registered trademark of the Hager & Meisinger GmbH, Germany

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# Benex®-Control Professional Benex®-Control Basic

Root Extraction System

85FLBM01 - 0315

# Benex®-Control Professional Benex®-Control Basic

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The Benex®-Control Bone Management® System is used for the controlled and safe extraction of tooth roots and tooth root fragments (e.g. after a trauma of the anterior teeth). Due to the innovative construction of the extractor, the root can be removed very easily and in an especially controlled manner without damaging any bone or soft tissue. This means that the optimal protection of the vestibular wall is ensured by means of the axial displacement of force.

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Separating the desmodontal fibers (Sharpey's fibers) in the sulcus using a periotome. Long and strong roots have to be axially loosened and luxated in both mesial and distal directions for a period of 30 seconds using a delicate lifting device / twister. In this process, care must be taken that no transversal movements are made. When extracting molars, the roots have to be separated before removal.

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
An axial drill hole should be made in the center of the tooth root and the dental root fragment that is still present using a suitable pilot bur of the appropriate size and a diamond-bit twist drill. This should have a depth of approx. 7 mm in hard tissue; a deeper drill hole is not necessary. In order to make twisting in the extraction screws easier, it is recommended to use a diamond-bit twist drill. Ideally, during the procedure, cooling should take place using water. In order to better remove the chips caused during drilling, the twist drill or pilot bur should be moved "intermittently".

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Depending on the indication, the right extraction screw can be carefully inserted into the tooth root using the right jollying aids. Note: Care must be taken that the jollying aids are pressed onto the extraction screw. This is the only way to ensure that there is a sufficient level of clamping force for safe transport.

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The images illustrate how the traction cord is threaded into the protection cap. Here it is of advantage to keep the traction cord very short.

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
The protection cap is now moved slightly upwards. Then the cord can be hooked in as described for the next step and the protective cap can be slid over the screw head.

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
The extractor is positioned on the adjacent tooth crowns with its tooth-sparing Teflon® surface. The opening of the round, rotating support should be in vestibular position here, whereby the extraction screw should be clearly visible. The traction cord is hooked into the extraction screw using the prepa-

red protection cap (see graphic), led through the pulley, and attached to the extraction slide hook. In the process, the screw and the cord must have the same axial direction.

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Once the extractor has been correctly positioned, extraction takes place by means of slow, controlled turning of the handscrew. It is important that the Parodontal fibers are tensioned with submaximal cable force at intervals of 30 seconds.

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
In order to guarantee an optimal vertical transmission of force, there must be a right-angled

plane present for application of the extractor. The quadrant support provides this right-angled application plane. Equally, a bad and uneven application surface can simply be corrected. The quadrant support makes it possible to evenly distribute the force applied in order to be gentle on teeth. Furthermore, it provides for an optimal application surface in the case of free-end situations and edentulous jaws.

- 1.) Fill the quadrant support with a rapidly hardening silicon or alginate and model a gap in the area of the buccal cavity.
- 2.) Position the quadrant support with the gap in the area of the root to be extracted at a right angle to the root axis.
- 3.) Remove the quadrant support when the silicon or alginate starts hardening. The entire hardening processing takes place extra-orally.
- 4.) Root luxation. Next, attach the Benex screw and reposition the quadrant support as well as the extractor.

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