# Inhalt Crestal-Lift-Control

Content Contenido

\* Crestal-Lift-Control Basic

	Vorkörner Initial bur Fresa inicial	Osteotomieinstrument Osteotomy instrument Fresa para osteotomia	Lagerfräser Ablative bur Fresa ablativa	Knochenmaterial-Verteiler Bonegraft condenser Distribuidores de material óseo	Pilotbohrer Pilot bur Fresa piloto
<b>CE 0044</b>					
<b>Fig.</b>	<b>186RF</b>	<b>859</b>	<b>TC21X</b>	<b>CL019</b>	<b>CL020 *</b>
Shank <sup>1</sup>	204	204	205	206	206
Size <sup>2</sup>	018	018	040	020	030
Length mm	12.0	10.0	14.0	2.5	2.5

<sup>1</sup> 204=RA, 205=RA L, 206=RA XL <sup>2</sup> Largest working part diameter in 1/10 mm

	Crestal-Bohrer Crestal drill Fresa crestal						Tiefenmess-Sonde Depth gauge Sonda de medición de profundidad
<b>Fig.</b>	<b>CL002</b>	<b>CL003</b>	<b>CL004 *</b>	<b>CL005 *</b>	<b>CL006 *</b>	<b>CL007 *</b>	<b>CL021</b>
Shank <sup>1</sup>	206	206	206	206	206	206	
Size <sup>2</sup>	028	031	033	036	038	041	
Length mm	2.0-12.0	2.0-12.0	2.0-12.0	2.0-12.0	2.0-12.0	2.0-12.0	illustrated 1:2 original total length 99.0 mm

<sup>1</sup> 206=RA XL <sup>2</sup> Largest working part diameter in 1/10 mm

	Stopphülsen Stop sleeves Casquillos de detención										
<b>Fig.</b>	<b>CL008</b>	<b>CL009</b>	<b>CL010</b>	<b>CL011</b>	<b>CL012</b>	<b>CL013</b>	<b>CL014</b>	<b>CL015 *</b>	<b>CL016 *</b>	<b>CL017 *</b>	<b>CL018 *</b>
Length mm	4.0	5.0	6.0	7.0	8.0	9.0	10.0	11.0	12.0	13.0	14.0
≙ max. drilling depth mm	12.0	11.0	10.0	9.0	8.0	7.0	6.0	5.0	4.0	3.0	2.0

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# Crestal-Lift-Control Crestal-Lift-Control Basic

85FL083 - 0115

Bone Management® is a registered trademark of the Hager & Meisinger GmbH, Germany

Hager & Meisinger GmbH | Hansemannstr. 10 | 41468 Neuss | Germany  
 Tel.: +49 (0) 21 31 20 12-0 | Fax: +49 (0) 21 31 20 12-222 | www.meisinger.de | info@meisinger.de  
 Meisinger USA, L.L.C. | 10200 E. Easter Avenue | Centennial, Colorado 80112 | USA  
 Tel.: +1 (303) 268-5400 | Toll free: +1 (866) 634-7464 | Fax: +1 (303) 268-5407  
 www.meisingerusa.com | info@meisingerusa.com

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
Crestal-Lift-Control is used for a minimally invasive, closed sinus-floor elevation within the scope of internal sinus lifts. Lifting the Schneiderian membrane already takes place during the transcrestal drilling procedure. The atraumatic design of the crestal drill serves to safely shape a conical bone flap. The head of the 4-blade crestal drill, which is curved inwards, can also be used to collect bone graft material in an ideal manner. In addition, the adapted stop sleeve system provides protection against membrane injuries.

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


The crestal drills are used in increasing diameters (up to the necessary diameter of the implant) for final drilling, applying the next largest stop sleeve (example here: CL014). This final drilling provides access to the maxillary sinus. Due to the atraumatic construction of the crestal drill with its 4 blades and the head which is inwardly curved, a conical bone flap can be shaped. This automatically causes the membrane to be raised by approx. 2 mm. Appropriate cooling has to be ensured during the entire drilling procedure. We recommend working intermittently at a max. of 800 min<sup>-1</sup>. In order to collect bone graft material, a speed of 50 min<sup>-1</sup> should not be exceeded.

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The mobility of the bone flap can be checked with the depth gauge. Combined with the next largest stop sleeve (example here: CL013), the membrane can be lifted by another millimeter if required.

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
Afterwards, the new sinus area, which had been formed by means of drilling, is filled with bone material.

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Using the bone graft condenser, the applied graft can be evenly distributed within the sinus cavity. We recommend using a bone graft condenser intermittently at a max of 10 min<sup>-1</sup>.

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- Residual bone height under 5 mm
- Not recommended in the case of complex morphology of the sinus floor (septum)
- Dentoalveolar growth which has not come to an end (exception: cases in which no dentoalveolar growth can be expected, e.g. ectodermal dysplasia)
- Active infections as well as pathological processes
- Diseases that impair bone metabolism

NOTE: A simple mucocele in the preparation area (sinus floor) does not represent any contraindication and can frequently be removed by means of aspiration.

CAUTION: Despite atraumatically designed instruments, care must be taken not to perforate the membrane.  
NOTE: Related contraindications and further instructions that are deemed set-specific and generally valid can be found at the download area of our website at „Application and Safety Instructions for Bone Management® Systeme“.

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