

Trephine Basic Kit

Trephine System for Extraction of Cylindrical Bone Grafts

Das Trephine Basic Kit bietet ein ausgewähltes Sortiment hochwertiger Trepane verschiedener Größen. Diese sind zum präzisen und einfachen Erstellen von zylindrischen Knochentransplantaten sowie zur Knochenentnahme geeignet. Die Trepane haben eine Tiefe von 10 mm und eignen sich speziell für den Einsatz in der Kieferchirurgie und Implantologie. Dank der gut erkennbaren Tiefenmarkierungen im Abstand von jeweils 2 mm kann beim Einsatz jederzeit leicht und sicher die Frästiefe abgelesen werden.

The Trephine Basic Kit provides a selected set of high-quality trephine drills in different sizes. These tools are designed specifically for creating cylindrical bone grafts quickly and accurately and for removing defined bone sections. The trephine drills have a depth of 10 mm and have been developed specifically with mandibular surgery and implantology in mind. Thanks to the well-defined laser depth marking in 2 mm intervals, the drilling depth can accurately be monitored.



Art.-No. 7120

