



Ost-Tray NEW

Selection of Manual Instruments for Oral Surgery | developed by Prof. Dr. Fouad Khoury

Das Bone Management® Ost-Tray bietet eine perfekte Zusammenstellung aus speziell konfigurierten manuellen chirurgischen Instrumenten. Neben einer durchdachten Auswahl an Handinstrumenten, wie hochwertigen Parodontalsonden, Küretten, Skalpellklingenhaltern, Mundspiegeln, großenabgestuften Mundkeilen, einer gebogenen Spülkanüle mit Olive sowie diversen Scheren, Nadelhaltern, Pinzetten, Klemmen, Zangen, Wundhaken und Zungenspateln, enthält das Kit insbesondere eine individuell konstruierte, kräftige Khoury-Sonde sowie auf Prof. Khoury abgestimmte Raspatorien und scharfe Löffel.

The Bone Management® Ost-Tray offers a perfect combination of specially configured manual surgical instruments. In addition to a well-thought-out selection of manual instruments, such as high-quality periodontal probes, curettes, scalpel blade holders, mouth mirrors, graduated mouth wedges, a curved irrigation cannula with bulb end and various scissors, needle holders, tweezers, clamps, forceps, retractors and tongue depressors, the kit contains an individually designed, strong Khoury-probe as well as raspatories and sharp spoons according to Prof. Khoury.



Art.-No. BOST0



Fig.	MI131	MI130
Name	Mundkeil / Mouth wedge	Mundkeil / Mouth wedge
Size	L	XL

Optional erhältlich:
Optionally available:

Spülkanüle mit Olive (Art.-No. MI164)
Irrigation Needle with Bulb End (Art.-No. MI164)

Mundkeil S (Art.-No. MI132)
Mouth wedge S (Art.-No. MI132)



Handinstrumente
Manual instruments



Fig.	MI056	MI057	MI070	MI071	MI077	MI139	MI156/165	MI160	MI163
Name	Scharfer Löffel 'Khoury', 2,0 mm Bone Curette 'Khoury', 2,0 mm	Scharfer Löffel 'Khoury', 3,0 mm Bone Curette 'Khoury', 3,0 mm	Raspatorium 'Molt' Periosteal Elevator 'Molt'	Raspatorium 'Khoury' Periosteal Elevator 'Khoury'	Skalpellklingenhalter Scalpel Handle	Kürette 'Columbia' Curette 'Columbia'	Mundspiegel Mouth Mirror	Parodontalsonde Periodontal Probe	Sonde, kräftig Probe, strong
Length mm	175.0	175.0	180.0	183.0	152.0	166.0	165.5	160.0	158.0

Handinstrumente
Manual instruments




Fig.	M1083	M1101	M1102	M1104	M1105	M1106
Name	Präparierschere 'Metzenbaum', TC Dissecting Scissors 'Metzenbaum', TC	Wundhaken 'Langenbeck' 42x10 mm Tissue Retractor 'Langenbeck' 42x10 mm	Wundhaken 'Langenbeck' 60x11 mm Tissue Retractor 'Langenbeck' 60x11 mm	Arterienklemme 'Halstead-Mosquito' Haemostatic Forceps 'Halstead-Mosquito'	Tuchklemme 'Tohoku' Towel Clip 'Tohoku'	Hohlmeißelzange 'Blumenthal' Bone Rongeur Forceps 'Blumenthal'
Length mm	180.0	220.0	220.0	145.0	105.0	155.0

Handinstrumente
Manual instruments



Fig.	M1041	M1048	M1050	M1158	M1089	M1098
Name	Nadelhalter 'Crile Wood', TC Needle Holder 'Crile Wood', TC	Pinzette Tweezers	Pinzette, TC Tweezers, TC	Pinzette 'College' Tweezers 'College'	Chirurgische Schere 'Kelly' Surgical Scissors 'Kelly'	Zungenspatel 'Brünnings' Tongue Depressor 'Brünnings'
Length mm	150.0	175.0	175.0	175.0	160.0	200.0